Run 1			
Static Data		Start Time (peltier on):	12:34 PM
Size of cold stage:	51 mm	End Time (peltier off):	13:29
Height of cold stage:	8 mm		
vacuum pressure (set)	50 Pa		
Z/TILT (distance of detector from stage)	5um		
vacuum current	12kv		
probe current	90		
Additional Comments	For this run I waited untill the v	accum pressure was 60P	a to start
Additional Comments	It pretty quickly turned to 50 Pa		
Kinetic Data			
Time	Action/observation	Temperature (actual)	Crystal Size
12:34	set temp to -35	25	
12:38	crystaks observed	-35	
12:41	set temp to -32 to take the pics		
12:48	took pics		
12:48	set to -25 to ablate	-26	
12:50	ablated to fast (set to -32)	-32.8	
12:54	set to -30 to see slower ablation		
12:56	slow ablation (pockmarks occu	-30.8	
	set -32 to take pics	-32.6	
	finished pics	-32.8	
1:02	set to -30 to see ablation	-30.8	
1:05	set to -31		
1:06	seems to be holding shape rela	-31.8	
	set to -33	-33.7	
	good steady re-growth	-33.7	
	set -34	-34.6	
	perfect re-growth	-34.7	
	growth roughening is appearing		
	set to -33 to slow for pic	01.0	
	set to -32 to slow even more		
	capturing immages of growth r	-32.7	
	set to -34	-34.6	
	ended experiment	-34.6	
1.20	опаса охроппон	01.0	
Run 2			
Static Data		Start Time (peltier on):	2:12 AN
Size of cold stage:	51 mm	End Time (peltier off):	2:34
Height of cold stage:	8 mm	Life Timo (poleor on).	2.0
vacuum pressure (set)	50 Pa		
Z/TILT (distance of detector from stage)	5um		
vacuum current	12kv		
probe current	90		
Additional Comments	vial was not compleaty frozen when I started		
Additional Comments	viai was not compleaty nozem	When I started	
Kinetic Data			
Time	Action/observation	Temperature (actual)	Crystal Size
		pressure is 140	Orysiai SIZE
	set temp to -35 temp is -27 and pressure is 80		reach)
	crystals observed (pressure 60		start
	crystals observed (pressure 60 crystal has good morphology		300x150um
			SOUX FOUUITI
	set to -33 to hold shape	-33.6	
	ant 22	-32.8	
2:21	set -32		
2:21 2:25	set -35 to try and grow new cry		
2:21 2:25 2:30		-34.7	