

Run 1				
Static Data			Start Time (peltier on):	1:01 AM
Size of cold stage:	51 mm		End Time (peltier off):	1:30
Height of cold stage:	8 mm			
vacuum pressure (set)	50 Pa			
Z/TILT (distance of detector from stage)	5um			
vacuum current	12kv			
probe current		85		
Additional Comments	I waited untill 60 Pa to trun on Peltier			
Additional Comments	The immagine seems different			
Kinetic Data				
Time	Action/observation		Temperature (actual)	Crystal Size
1:01	set temp to -37		25 C	
1:04		-26		
1:05	crystals observed		-36.6	
1:06	vacuum pressure went down to 40 Pa instead of 50			
1:08	temp set to -35		-37.7	
1:11	Took 3D immagine (case 1.0)		-35.8	
1:13	set temp to -30 for ablation		-35.6	
1:14	Took 3D immagine (case 1.1)		-30.9	
1:16	Took 3D immagine (case 1.2)		-30.9	
Note: Fast er ablation seems to be occuring and this is causing small pockmarks				
1:18	set temp to -32 to try and slow ablation			
Note: crystal seems to be visibly dissappearing				
1:19	Took 3D immagine (case 1.3)			
1:21	I zoomed in on one of the prismatic facets and took a 3D immagine (case 1.4)		-32.9	
1:23	I zoomed in on one of the same prismatic facet but farther down and took at 3D immagine (case 1.5)		-32.9	
1:27	Took 3D immagine of whole ice crystal (case 1.6)		-32.9	
1:28	Ended experiment		-32.8	
Run 2 (this is using the normal SEM and X-Ray crystaloraphy of filter)				
Static Data				
Size stage:	15mm			
Height of cold stage:		-5		
vacuum pressure (set)	40 Pa			
Z/TILT (distance of detector from stage)	10mm			
vacuum current	17kV			
probe current		50		
Additional Comments	I was working on figuring out all of the settings so did not write as much down			
Additional Comments				
Kinetic Data				
Time	Action/observation		Temperature (actual)	Crystal Size