Run 1 (K2SO4)				
Static Data		Start Time (peltier on):	12:13 PM	
Size of cold stage:	51 mm	End Time (peltier off):		
Height of cold stage:	8 mm			
vacuum pressure (set)	50 Pa			
Z/TILT (distance of detector from stage)	5um			
vacuum current	12kv			
probe current	85			
Additional Comments	pressure is 50 Pa			
Additional Comments				
Kinetic Data				
Time	Action/observation	Temperature (actual)	Crystal Size	
2:4	temp set to -34	25		

Note: very little data available here due to the fact that more attention was being paid to taking pictures and in person observations (this was the first run of this done)